

SaeboStim Go



Close the Gap. Wireless wearable stimulation for shoulder subluxation.

The **SaeboStim Go**, a wireless electrical stimulation device, is designed to prevent or reduce subluxation by strengthening the weakened muscles responsible for the instability or dislocation. The affordable and compact device works by delivering electrical signals through its wireless electrodes to the targeted shoulder muscles that lift the arm back into the joint.

The **SaeboStim Go** complements the Saebo program. Whether applying the stimulation before or after the Saebo treatment or combining stimulation with other Saebo technology, providing multiple evidence-based treatment solutions can enhance overall functional recovery.

Features:

- Affordable and simple to use.
- Wireless, sleek, and lightweight.
- Reusable gel pads are made of the highest quality gel available, and due to the design of the device, are a low cost when you need to refill.
- Includes only 1 preset program, so no guess work.
- Rechargeable battery.
- Utilizes a patent-pending ultra-thin, rotating wing structure that is fully flexible which allows you to apply device to specific areas.
- Laser printed silver conductive surface to provide an unmatched quality of stimulation that doesn't feel sharp like other devices.

Indications of Use:

- Shoulder subluxation
- Disuse atrophy/weakness
- Muscle re-education
- Range of motion

Not just for Subluxation!

The **SaeboStim Go** can be used for other parts of the body as well. The device's program parameters are clinically designed to improve function, slow down atrophy, and strengthen weakened muscles in both upper and lower extremities.

